The Glutavac test

GB

What is Clutavac?

Glutavac is a trade name registered for evacuated test tubes containing a 1.4% glutaraldehyde solution.

The Glutavac test

The Glutavac test is a glutaraldehyde test suitable for use on whole blood for the diagnosis of acute inflammation (increased plasma fibrinogen concentration) and chronic inflammatory infections (increased plasma immunoglobulin concentration) in large animals (cattle, small ruminants, horses and pigs).

Principle

Glutaraldehyde reacts rapidly with protein amino groups, and cross-links the protein molecules irreversibly. At adjusted concentrations of glutaraldehyde, there is an inverse relationship between fibrinogen plus immunoglobulin concentrations, and the time taken for cross-linking leading to gel formation. In normal cattle, pigs, and horses, gel formation takes more than 15 minutes.

Note

Not suitable for patients with severe anaemia or severe dehydration. Not suitable for animals younger than 6 months of age. Test tubes which do not draw blood to the red marks on the label should be discarded.

Experiences

The glutaraldehyde test on whole blood was first developed by Sandholm (1974, 1974a) as a rapid cow-side test for the diagnosis of chronic infections, such as those involved in traumatic peritonitis. Later investigations showed it equally useful in acute infections with increased fibrinogen concentration (Liberg, 1978). This was confirmed in experimentally induced focal ruminitis by Jørgensen, Basse, and Aslan (1990). Nielsen (1975) found the test valuable in both cattle and horses, and Hansen (1985) applied it in pigs, and showed that emergency slaughtered cows were accepted at meat inspections when the test was negative, whereas 62% of cows with a positive test were condemned. In his thesis, Martens (1977) examined 79 healthy and 112 sick cattle. He concluded that the glutaraldehyde test proved to be especially easy and rapid and that it is also suitable for routine use in the veterinary surgery. Liberg (1981) examined glutaraldehyde and formal-gel tests in 82 cattle with traumatic peritonitis. According to their reaction to the tests, the cases were divided into 6 reaction groups which were considered to reflect different intensities and/or stages of the disease. In general, large animals with acute and chronic inflammation will give a positive test (gel formation) within 15 minutes. The test will remain negative in cases of displaced organs, superficial gastro-intestinal infections and other conditions with no or low inflammatory response.

The Glutavac test gives you the answer within 15 minutes.

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