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
Left blank intentionally

6. Specifications

Weight	: approx. 2 kg
Storage temperature	: -20°C to +50°C
Operating temperature	: 0-40°C
Humidity	: 0-100%
Pressure	: 630-1060 mbar
Classification according to EU Medical Device Directive 93/42/EEC	: IIb
Mode of operation	: Continuous


Left blank intentionally


Description	Part no.	Qty.
1. Absorber valve	11900-120	1
2. Valve dome with O-ring	11890-40	2
3. Valve flap with seat	11890-31	2
4. DISC	12360-06	1
5. Silicon membrane	12360-08	1
6. Top w/automatic	11902-30	1
7. O-ring 22 mm	36267	6
8. O-ring for absorber valve	36289	1
9. Gasket for absorber valve	36511	1
10. Spring	11890-13	1
11. Water trap	11900-205	1
12. Top valve	11900-112	1
13. Locking	36737	2
14. Spacer	11900-111	1
15. Silicone O-ring	36258	1

Service kit for CUBE III:
 With 45 mm cones P/N 36428
 With 50 mm cones P/N 36429

- Part list:**
- 1. 45 mm cone11054-04
 - 2. O-ring EPDM36274
 - 3. Chamber11054-03
 - 4. Gasket11053-07
 - 5. Lower part11053-02B
 - 6. O-ring EPDM36303

5.10 Service Parts CUBE

5.11 Part List CO₂ Absorber

5.7 Cleaning Methods

5.7.1 Washing machinery

The washing machinery used must be constructed to allow for the washing of instruments and anaesthesia equipment. The washing programme and detergents used must be suitable for washing rubber and plastic parts max. Ph value 10. The washing procedure should be finalized by rinsing in clean water, heat disinfection and drying.

Warning:

When removing the CUBE from washing machine some hot water might be into the channels, because of scalding.

5.7.2 Chemical disinfection

Observe the manufacturer's instructions carefully. Disinfectants must be suitable for the components for which they are used. Do not expose components to disinfectants for longer time than needed. Rinse all components in sterile water after disinfection and dry at room temperature.

Note:

Some disinfectants, such as phenol and aldehyde compounds, can be absorbed and accumulated in plastic and rubber parts, when disinfection has taken place repeatedly.

Note:

The lifetime of rubber parts can be extended considerably by avoiding at decontamination and sterilization at high temperatures.

5.8 Safety Control

- After cleaning and sterilisation, reassemble and verify proper function of the system.
- Do not add components to the system without making sure they are clean and in good condition.
- Do not use the system before a leakage- and function test has been carried out.

5.9 Maintenance

Preventive maintenance shall be carried out at least once a year by a qualified technician.

1. Change the 22 mm O-rings on the INSP, EXP, VENT and ABSORBER RETURN cones and valve seats.

2. Change valve flaps incl. valve seats and the valve domes in the in/exp. Valves.

3. Change the gasket and O-ring in the absorber valve by removing the screw in the bottom of the CUBE, removing the absorber valve with spring and disassemble the valve by removing the lock-ring.

4. Change the 22 mm and 30 mm O-rings on the respiration valve.

5. Change the silicon membrane in the change-over valve on type 11902/11907.

6. Change O-rings and gasket on the absorber.

7. Change silicone O-ring on watertap.

8. Disassemble change-over valve on type 11900/11905 and 11902/11907 (use special tool P/N 72210-10). Clean the valve and change O-rings.

5.4.1 Assembly of the CUBE Circle

Before assembly all parts must be completely dry, and the O-rings on the absorber valve and INSP/EXP valves should be coated with a thin film of special grease Dameca P/N 36825.

- Mount the insp/exp valves (5) with the groove in the valve body fitting the pin in the CUBE. Check that the valves are correct mounted by carefully pushing the valves in. Check that the membranes have no visible damages. Change the membranes if they are damaged (use valve P/N 11890-31). Don't touch the valve flap.
- Mount the valve caps (4).
- Mount the absorber valve (3) with spring (2). On type 11902/11907 put the metal plate (6) into the valve housing, mount the silicon membrane (7), and screw the lid (8) on the housing.

After washing a thin film of special grease (P/N36825) shall be applied to the absorber valve and O-rings, except the upper O-ring in the absorber cone. Otherwise the valve can be stuck, and O-rings crack.

Note:

Special grease (P/N36825) can expose dangerous vapour, if heated to 300-350°C, avoid smoking when using special grease

Warning:

The circle can now be mounted as described in section 4.1.

5.5 CO₂ Absorber


Fig. 5-2 CO₂ absorber

Disassemble the components of the absorber as shown in the figure, and clean and sterilize as other patient related parts. Let all parts dry thoroughly before assembly and apply a thin coating of special grease (DAMECA P/N 36825) to the O-rings in order to ensure good connection.

Note:

Before reassembling visually check that all components are in correct condition and completely dry.

5.6 Hoses, Y-piece and Manual Breathing Bag

Clean, disinfect and dry in accordance with normal hospital routines.

5.4 CUBE Circle System


Fig. 5-1 Cube circle type 11902 with pneumatic change-over valve

Before maintenance and cleaning, disassemble the parts of the CUBE circle system as follows:

- Remove the screw/water trap (1), and remove the spring (2) with absorber valve (3).
- Remove the caps (4) for the insp/exp valves (5), and remove the valves.
- Remove the cover (8) of the change-over valve.
- Remove the silicon membrane (7), and the metal disc (6).

Note:

The change-over valve in type 11900 shall not be taken apart for cleaning. Type 11901/11904 is not equipped with a change-over valve.

The CUBE circle and the parts can now be cleaned and sterilized.

5. Cleaning and Maintenance

5.1 Routine Cleaning

Parts of the patient breathing system can be washed in a washing machine.

5.2 After An Infected Case

Soak all parts of the patient breathing system in a disinfectant.

Then wash according to the instructions for routine cleaning

5.3 After Cleaning

Check visually that all parts are totally dry undamaged and connected properly.

Note: Some disinfectants, such as phenol and aldehyde compounds, can be absorbed and accumulated in plastic and rubber parts, when disinfection has taken place repeatedly

Note: When cleaning in washing machine the detergent used must comply with the equipment. The detergent shall have Ph values between 7-10 in final solution.

Note: All parts should be completely dry before reassembling. Apply a thin coating of special grease to the O-rings of the absorber valve and INSP/EXP valves use Dameca part nr 36825.

Warning: Do not lubricate any part with any type of combustible lubricant as they may cause a danger of explosion. Use special grease (Dameca part nr 36825) only.

Note: Before the parts are re-assembled they shall be completely dry and in good condition.

Warning: If the system is used for more than one patient without cleaning, precautions to avoid cross infection must be taken.

Warning: When handling the Cube body after cleaning, beware of residual hot water in the internal channels.

When the system is assembled and before connecting it to a patient it must be checked for integrity and function.

4.2 Integrity Check (tightness)

1. Open the vaporizer.

2. Close the regulators on the gas box.

3. Close the respiration valve in the patient breathing system.

4. Place the patient connector (Y-piece or similar) on the test gauge.

5. Open the O₂ regulator, until the pressure reading on the test gauge is 30 cm H₂O.

6. Adjust the O₂ regulator, until the pressure reading on the test gauge is stable at 30 cm H₂O.

7. Read the O₂ rotameter.

8. If the flow is higher than 150 ml/min there is a leak. Locate the leak, repair it, and repeat the test.

Note:

The absorber shall be connected and open during the test.

4.3 Function Test

1. Open the spill valve (Bernier valve setting = 20cmH₂O), and block the patient connection.
2. Adjust the O₂ flow to 4L/min. and fill the system.
3. Squeeze the respiration bag a couple of times and check that the spill valve opens, and the inspiration/expirator valves are working.

Note:

Standard microbiological filters, gas sample adapters and adapters for flow measurement can be used with the system.

Rubber articles should frequently be checked for cracks or holes, as these are often causing leaks.


Note:

The maximum allowable leakage for the dose system is 50 ml/min., and for the patient breathing system 100 ml/min. At a pressure of 30cmH₂O.

By following the procedure described, the leakage in both dose- and patient-system is tested.

4. Operation

4.1 Assembly

qnsfob29;p20061.wmf

Note: It is recommended not to use heated humidifiers with rebreathing patient breathing circuits. If a heated humidifier for any reason is used, it must be checked frequently that it does not overheat.

Warning: The use of antistatic or electrically conductive breathing tubes when utilizing high frequency surgery equipment may cause burns and is therefore not recommended in any application of this machine.

Warning: The system must not be changed, or modified with components not specified by Dameca.

3. Function


Fig. 3-1 Circle system with absorber valve


Fig. 3-2 Circle system without absorber valve

The CUBE circle system is a semi-closed patient breathing system with one-way valves in both the inspiratory and the expiratory limb, and an absorber for eliminating CO₂.

The system is suitable for high flow (approx. the patients minute volume), medium flow (approx. 1 - 1.5 ltr/min) and low flow (approx. 0.5 - 1 ltr/min) of freshgas for patient ventilation.

The system is available with or without an absorber valve. If an absorber valve is present, it allows changing the absorber during use, without letting anaesthetic gases into the operating room.

Changing absorber during an operation:

1. Remove the absorber from the Cube circle.
- Note:** Leave the top cone and tubing on the absorber.
2. Remove the top cone and tubing from the old absorber, and place it on the new absorber.
 3. Place the new absorber (with top cone and tubing) in the cone on the Cube circle, by pushing it down and twisting it slightly.

2.3 Failure of Operation

- If the equipment fails to respond as described, do not use until the situation has been corrected by qualified service personnel.

2.4 Occupational Hazard

- Do not use this equipment in insufficiently ventilated areas. Use scavenging system to avoid increased concentrations of anaesthetic gases in operating rooms.

2.5 Cautions

- This equipment should be serviced at least once a year by authorised service personnel.
- Use only original **DAMECA** sampling lines and accessories; other sampling lines may cause inaccurate readings and malfunctions.

2.6 Equipment Safety Symbols

Attention consult the accompanying documents.

Note:
Tubings used with this product must fulfill EN 12342 standard, and respirations bags must fulfill EN 1820 standard.

- Use only original **DAMECA** O-rings; other O-rings may cause problems if they are not of the same size or material.

2. Safety Precautions

Read this manual carefully and familiarize yourself with the machine before you start to operate the machine.

Always check life-supporting equipment for proper function before use.

A WARNING indicates a situation in which the user or the patient may be in danger of injury death.

A CAUTION indicates a situation that may lead to equipment damage or malfunction.

A NOTE indicates a point of particular interest for more efficient and convenient operation.

WARNINGS:

2.1 Explosion Hazard

- Do not use this machine in the presence of flammable anaesthetics. Only anaesthetics classified according to IEC 60601-2-13 as non flammable must be used.

- No oil or grease should be used in this machine. This also applies to silicone based lubricants and inflammable oil or grease. Only PTFE based oxygen compatible lubricants must be used, (e.g. DAMECA special grease P/N 36825).

2.2 Patient Safety

- It is required that anaesthesia workstations shall always be used with monitors for O₂, anaesthetic agent, CO₂ concentrations, and expired volume.

- Do not use the equipment for a clinical procedure without first ensuring that the complete system is operating correctly.

- Ensure that all connections are secure and check for proper operation of the machine before use. In particular the breathing system, reusable adapters, and O-rings should be checked to ensure that connections are not loose and the breathing circuit is not leaking.

- The use of antistatic or electrically conductive breathing tubes when utilizing high frequency surgery equipment may cause burns and is therefore not recommended in any application of this machine.

- Disposal airway adapters and sampling tubes are for single use only. Re-useable adapters should be sterilized according to the hospitals normal routines.

- Constant attention by a qualified individual is required whenever a patient is under anaesthesia or connected to a ventilator. Some malfunction of equipment require immediate action. A malfunction may pass unnoted in spite of equipment alarms.

- The change of colour of the soda lime depend of the different manufacturers. See the manufacturers instructions.

- Patient breathing tubings used with Cube III must comply with EN 3993-1 or similar standard.
- Gas sample tubings used with Cube III must comply with EN 1707 or similar standard.

Left blank intentionally

1. Introduction

Dameca CUBE circle is a rebreathing patient

breathing system.

The system consist of an integrated block with inspiration- and expiration valves, BAG/VENT

change-over valve, holder for the system,

transparent CO₂ absorber, manual ventilation

bag, mask and tubings.

When changing the absorber an automatic valve

will close in order not to loose any fresh gas.

Can be ordered without this valve.

All parts can be washed in washing machine.

CUBE III is made from plastic material (POM).

1.1 Versions

There are four different systems with 45 or 50 mm absorber cone:

System	Part number
CUBE circuit with manual change-over valve, 45 mm/50 mm	11900/11905
CUBE circuit with pneumatic change-over valve, 45 mm/50 mm	11902/11907
CUBE circuit without change-over valve, only for VENT, 45 mm/50 mm	11904/11909
CUBE circuit without change-over valve, only for BAG, 45 mm/50 mm	11901/11906

1.2 Intended Use

The patient breathing system is intended to be used with Dameca anaesthesia machines or machines that fulfill EN740 standard or similar standards.

The system is intended for use with adult and paediatric patients, but not to be used with neonatal patients.

The system must be used by qualified anaesthesia personal only.

1.3 Performance

Expiratory and inspiratory resistance at a flow of:

60 ltr/min	30 ltr/min	5 ltr/min
Expiratory resistance (x100 Pa)	4.5	2.0
Inspiratory resistance (x100 Pa)	4.5	2.0
Expiratory resistance (x100 Pa)	0.3	0.3
Inspiratory resistance (x100 Pa)	0.3	0.3

Compliance = 3 ml/hPa with adult tubings (ø22mm)

Note:

The system cannot be used with Dameca re-useable absorbers produced earlier than medium 1999. After 1999 the re-useable absorber male cone is changed to an inner diameter of 32 mm. Early produced CO₂ absorbers can be modified by changing the lower part: P/N 11054-10 with 45 mm cone and P/N 11053-10 with 50 mm cone.

Left blank intentionally

Contents

1. Introduction.....	1-1
1.1 Versions.....	1-1
1.2 Intended Use.....	1-1
1.3 Performance.....	1-1
2. Safety Precautions.....	2-1
2.1 Explosion Hazard.....	2-1
2.2 Patient Safety.....	2-1
2.3 Failure of Operation.....	2-2
2.4 Occupational Hazard.....	2-2
2.5 Cautions.....	2-2
2.6 Equipment Safety Symbols.....	2-2
3. Function.....	3-1
4. Operation.....	4-1
4.1 Assembly.....	4-1
4.2 Integrity Check (tightness).....	4-2
4.3 Function Test.....	4-2
5. Cleaning and Maintenance.....	5-1
5.1 Routine Cleaning.....	5-1
5.2 After An Infected Case.....	5-1
5.3 After Cleaning.....	5-1
5.4 CUBE Circle System.....	5-2
5.5 CO ₂ Absorber.....	5-3
5.6 Hoses, Y-piece and Manual Breathing Bag.....	5-3
5.7 Cleaning Methods.....	5-4
5.8 Safety Control.....	5-4
5.9 Maintenance.....	5-4
5.10 Service Parts CUBE.....	5-5
5.11 Parts List CO ₂ Absorber.....	5-5
6. Specifications.....	6-1

Left blank intentionally

User Manual for Dameca CUBE III Circle Model 11900

Important !

We are constantly developing our products, and consequently reserve the right to make alterations in design and equipment without prior notice.

DAMECA a/s issue a guarantee against manufacturing and material defects for a period of one year from the date of delivery. Excluded from this guarantee are rubber articles and anaesthesia valves with plastic parts.

DAMECA a/s will repair or replace all defective components during the period of guarantee, provided (1) that the unit has been used correctly in accordance with the instructions for use and (2) that the instructions for maintenance have been observed carefully.

Defective products or components must be returned to **DAMECA a/s** with a detailed description of the defect in question.

If repairs are carried out by anyone other than **DAMECA a/s** or by a **DAMECA a/s** appointed agent, the guarantee no longer applies.

Nobody but **DAMECA a/s** can issue any other form for guarantee, either in writing or orally.

DAMECA a/s has no responsibility for any translation of this English version of the manual.


Preventive maintenance of this unit shall be carried out at least once a year either by **DAMECA a/s** or an authorized technician.

Document No. 11900-90
August, 2004

Conformity according to the
Council Directive 93/42/EEC
concerning Medical Devices.

Dameca CUBE III Patient System Model 11900

Dameca