ASP™

REF 20391

CIDEX™ OPA

- en** CIDEX OPA - *ortho*-Phthalaldehyde High Level Disinfectant Solution
- de** CIDEX OPA - Hochwirksame Desinfektionslösung mit dem Wirkstoff *ortho*-Phthalaldehyd
- fr** CIDEX OPA - Solution désinfectante d'*ortho*-phtalaldéhyde
- it** CIDEX OPA - Soluzione disinfettante ad alto livello a base di *orto*-ftalaldeide
- es** CIDEX OPA - *orto*-ftalaldeído, Solución de Gran Poder Desinfectante
- sv** CIDEX OPA - *orto*-ftalaldehydlösning, höggradigt desinfektionsmedel
- nl** CIDEX OPA - *ortho*-ftaalaldehyde-oplossing high level-desinfectiemiddel
- pt** CIDEX OPA - *orto*-ftalaldeído Solução Desinfectante de Alto Nível
- fi** CIDEX OPA - *orto*-ftalaldehydiä sisältävä desinfiointiliuos korkeatasoiseen desinfiointiin
- da** CIDEX OPA - *ortho*-phthalaldehyd højniveau desinfektions opløsning
- no** CIDEX OPA - *ortho*-ftalaldehyd høyvirksom desinfeksjonsvæske (HLD)
- el** CIDEX OPA - ορθο-φθαλάλδεϋδι Αδσθεομαίθευ Αέϋθεοία Υπερειύ Αθεοϋαίθε
- pl** CIDEX OPA - Roztwór silnie dezynfekujący na bazie aldehydu *orto*-ftalowego
- hu** CIDEX OPA - *Orto*-ftálaldehyd magas szintű fertőtlenítő oldat
- cs** CIDEX OPA - Vysoce účinný dezinfekční roztok s obsahem *ortho*-ftalaldehydu
- sk** CIDEX OPA - *orto*-ftalátaldehyd roztok dezinfekčného prostriedku s vyšším dezinfekčným účinkom
- et** CIDEX OPA - *orto*-ftaalaldehyüüdi tugevalt desinfitseerivate omadustega lahus
- tr** CIDEX OPA - *ortofitalaldehyt* Yüksek Seviye Dezenfeksiyon Solüsyonu
- ru** CIDEX OPA - Раствор орто-фталальдегида для дезинфекции высокого уровня

English	4
Deutsch	13
Français	23
Italiano	33
Español	43
Svenska	53
Nederlands	62
Português	72
Suomi	82
Dansk	92
Norsk	101
Ελληνικά	110
Polski	121
Magyar	131
Čeština	141
Slovenčina	151
Eesti	161
Türkçe	170
Русский	179

INSTRUCTIONS FOR USE

CIDEX™ OPA *ortho*-Phthalaldehyde Solution High Level Disinfectant

Active Ingredient <i>ortho</i> -phthalaldehyde	0.55%
Inert Ingredients	99.45%
Dipotassium hydrogen phosphate	
Potassium dihydrogen phosphate	
Benzotriazole	
Citric acid	
D&C Green Dye #5	
N-(hydroxyethyl)-ethylenediaminetriacetic acid (HEDTA)	
Water	
Total	100.00%

Does not require activation before use.

INDICATIONS FOR USE

CIDEX OPA Solution is a high level disinfectant for reprocessing heat sensitive medical devices used in the medical area, when used according to the Directions for Use. CIDEX OPA Solution is intended for use in manual (bucket and tray) systems made from polypropylene, acrylonitrile-butadiene-styrene (ABS), polyethylene, glass-filled polypropylene and/or polycarbonate plastics. CIDEX OPA Solution may also be used in automated endoscope reprocessors according to the manufacturer's instructions and should be monitored with CIDEX OPA Solution Test Strips.

See DIRECTIONS FOR USE – Reusage for Disinfection.

Medical devices reprocessed in CIDEX OPA Solution must first be cleaned according to a validated cleaning protocol or standard.

Antimicrobial Level of Activity: CIDEX OPA Solution can be used at the following antimicrobial level of activity:

High Level Disinfectant: CIDEX OPA Solution is a high level disinfectant for semi-critical devices when used or reused, according to the Directions for Use, at or above its Minimum Effective Concentration (MEC) of 0.3% as determined by CIDEX OPA Solution Test Strips, at 20°C (68°F) with an immersion time of at least 5 minutes for a reuse period of up to 14 days.

Reuse Period for Disinfection: CIDEX OPA Solution has demonstrated disinfection efficacy in the presence of 5% organic soil contamination and microbiological burden during reuse. CIDEX OPA Solution may be reused for up to a maximum of 14 days provided the required conditions of *ortho*-phthalaldehyde concentration and temperature exist based upon monitoring described in the Directions for Use. DO NOT rely solely on days in use. Concentration of this product during its reuse life must be verified by the CIDEX OPA Solution Test Strip prior to each use to determine that the concentration of *ortho*-phthalaldehyde is above the MEC of 0.3%. The product must be discarded after 14 days, even if the CIDEX OPA Solution Test Strip indicates a concentration above the MEC.

General Information on Selection and Use of Disinfectants for Medical Device Reprocessing: Choose a disinfectant with the level of antimicrobial activity that is appropriate for the reusable device. Follow the reusable device labeling and standard institutional practices. In the absence of complete instructions, use the following process:

First, for patient contacting devices, determine whether the reusable device to be reprocessed is a critical or semi-critical device.

Critical device: *Presents a high risk of infection if not sterile. Routinely penetrates the skin or mucous membranes during use or is otherwise used in normally sterile tissue of the body.*

Semi-critical device: *Makes contact with mucous membranes but does not ordinarily penetrate normally sterile areas of the body.*

Second, determine if sterilization or high level disinfection is required.

Critical device (e.g. laparoscopes and microsurgical instruments): *Sterilization is required.*

Semi-critical reusable device (e.g. gastroscopes): *Sterilization is required whenever feasible; where not feasible, high level disinfection is the minimum acceptable process.*

Third, select a disinfectant that is labeled for the appropriate antimicrobial level and is compatible with the reusable device. Follow directions for the disinfectant.

CIDEX OPA was tested following disinfection by immersion as stated per the directions for use at 5 minutes, 20°C with MEC of 0.3% under clean conditions, where test results meet all requirements (per EN 14885:2018) for the Spectrum of Activity listed in the table below.

<u>Type of Activity</u>	<u>Applicable Standard</u>
Bactericidal	EN 13727, EN 14561
Yeasticidal	EN 13624, EN 14562
Fungicidal	EN 13624, EN 14562
Tuberculocidal	EN 14348, EN 14563
Mycobactericidal	EN 14348, EN 14563
Virucidal	EN 14476

The following list includes additional microorganisms tested for CIDEX OPA Solution:

FUNGI

Trichophyton mentagrophytes

VIRUSES NON-ENVELOPED

Poliovirus Type 1
Rhinovirus Type 42
Adenovirus Type 2
Coxsackievirus Type B-3

VIRUSES ENVELOPED

Coronavirus
Cytomegalovirus
Influenza Virus [Hong Kong]
HIV-1
Herpes simplex Types 1,2
Vaccinia (Wyeth)

MICROORGANISM VEGETATIVE MICROORGANISMS

Staphylococcus aureus
Salmonella enterica (formerly
Salmonella choleraesuis)
Pseudomonas aeruginosa
Mycobacterium bovis

Material Compatibility: CIDEX OPA Solution has been tested and found to be compatible with the materials shown below.

METALS¹

Aluminum
Anodized aluminum²
Brass
Carbon steel
Chrome plated brass²
Chrome plated steel²
Copper
Nickel plated brass²
Nickel silver alloy²
Stainless steel³
Titanium
Tungsten carbide²
Vanadium steel⁴

PLASTICS⁵

Polymethylmethacrylate (Acrylic)
Nylon
Polyethylene terephthalate (Polyester)
Polystyrene
Polyvinylchloride (PVC)⁶
Acrylonitrile/butadiene/styrene (ABS)
Polysulfone
Polycarbonate⁷
Polyethylene
Polypropylene
Acetal
PTFE
Polyamide

ELASTOMERS⁵

Polychloroprene (Neoprene)
Kraton G
Polyurethane
Silicone rubber⁶
Natural rubber latex

ADHESIVES⁵

Cyanoacrylate⁸
EPO-TEK 301 epoxy⁸
EPO-TEK 353 epoxy

DENTAL MATERIALS⁹

Polysulphide
Addition silicone
Polyether

1. Exposed to 31 days (744 hours) of continuous contact with CIDEX OPA Solution with no effect unless otherwise noted.
2. Shows signs of surface discoloration at 7 days or greater.
3. Most grades tested show no effect. Others may exhibit slight discoloration at 7 days or greater. Stainless steel 440 shows rust at 14 days immersion.
4. Treated with 500 cycles of CIDEX OPA Solution. Surface breakdown noted after 150 cycles (25 hour total contact).

5. Exposed to 7 days of continuous contact with CIDEX OPA Solution with no effect unless otherwise noted.
6. Some grades or applications exhibit discoloration.
7. Some sonic welded parts may exhibit crazing.
8. Some loss in shear strength but show no signs of severe degradation.
9. Dental materials treated with 1 cycle of CIDEX OPA Solution for 10 minutes with no effects.

Fujinon, Olympus and Pentax endoscopes are compatible with CIDEX OPA Solution. If questions arise regarding the compatibility of a device with CIDEX OPA Solution, contact the device manufacturer.

Cleaning Agent Compatibility: CIDEX OPA Solution is compatible with enzymatic detergents which are near neutral (6-8) in pH, low foaming, and easily rinsed from equipment (e.g., CIDEZYME™ and CIDEZYME™ LF Enzymatic Detergents). Detergents that are either highly acidic or alkaline are not recommended as cleaning agents.

CONTRAINDICATIONS

1. CIDEX OPA Solution should not be utilized to process any urological instrumentation used to treat patients with a history of bladder cancer. In rare instances CIDEX OPA Solution has been associated with anaphylaxis-like reactions in bladder cancer patients under going repeated cystoscopies.
2. CIDEX OPA Solution should not be utilized to process instrumentation for patients with known sensitivity to CIDEX OPA Solution or any of its components.
3. CIDEX OPA Solution should not be used to sterilize heat sensitive medical devices.

WARNINGS

CAUTION:

Contains 0.55% *ortho*-Phthalaldehyde.
May produce an allergic reaction, including anaphylaxis.

Hazard statement

May cause an allergic skin reaction. Toxic to aquatic life with long lasting effects.

Precautionary statement

Avoid breathing mist/vapours. Wear protective gloves. Contaminated work clothing should not be allowed out of the workplace. Avoid release to the environment.
IF ON SKIN: Wash with plenty of water. If skin irritation or rash occurs: Get medical advice/attention. Take off contaminated clothing and wash it before reuse. Collect spillage. Store away from incompatible materials. Dispose of contents/container in accordance with local/regional/national/international regulations.
Contains Phthalaldehyde. May produce an allergic reaction.

1. May elicit an allergic reaction. Possible allergic reactions, including anaphylaxis have been reported in rare instances. In the majority of these instances health care workers were not using the product in a well-ventilated room or not wearing proper personal protective equipment. (See PRECAUTIONS).
2. Avoid contact with eyes, skin, or clothing. (See PRECAUTIONS – for important information on how to protect eyes, skin and clothing.) Direct contact with eyes may cause irritation. Direct contact with skin may cause temporary staining. Repeated contact may cause skin and respiratory sensitization. Sensitization can result in anaphylaxis in rare instances. In case of eye contact, immediately flush eyes with large quantities of water for at least 15 minutes. Seek medical attention. In case of skin contact, immediately wash with soap and water. Refer to the SDS for additional information. Do not form sprays, mists or aerosols of this product.
3. Ingestion may cause irritation or chemical burns of the mouth, throat, esophagus and stomach. If swallowed, DO NOT INDUCE VOMITING. Drink large quantities of water and call a physician immediately. Probable mucosal damage from oral exposure may contraindicate the use of gastric lavage.
4. Avoid exposure to *ortho*-Phthalaldehyde vapors, as they may be irritating to the respiratory tract and eyes. May cause stinging sensation in the nose and throat, discharge, coughing, chest discomfort and tightness, difficulty with breathing, wheezing, tightening of throat, urticaria (hives), and/or sensitization/anaphylaxis in rare instances, rash, loss of smell, tingling of mouth or lips, dry mouth or headache. May aggravate a pre-existing asthma or bronchitis condition. In case of adverse reactions from inhalation of vapor, move to fresh air. If breathing is difficult, oxygen may be given by qualified personnel. If symptoms persist, seek medical attention.
5. The use of CIDEX OPA Solution with semi-critical devices must be part of a validated rinsing procedure as provided by the device manufacturer. See DIRECTIONS FOR USE Rinsing Instructions for important information on rinsing. Improper rinsing can result in anaphylaxis in rare instances.
6. ALWAYS follow the Directions For Use Rinsing Instructions (Part B) and the SPECIAL INSTRUCTIONS for transesophageal echocardiography (TEE) probes in Part C EXACTLY or residues of CIDEX OPA may remain on the device. Failure to follow rinsing instructions exactly has resulted in reports of chemical burns, irritation, and staining of the mouth, throat, oesophagus and stomach. Improper rinsing can result in anaphylaxis in rare instances.

PRECAUTIONS

Follow Hospital policy and protocol when handling and cleaning soiled devices.

1. When disinfecting devices, use gloves of appropriate type and length. Nitrile rubber gloves or butyl rubber gloves may be used. When using latex rubber gloves, the user should double glove and/or change single gloves frequently, e.g., after 10 minutes of exposure. For those individuals who are sensitive to latex or other components in latex gloves, 100% synthetic copolymer gloves are recommended.
Note: Contact with CIDEX OPA Solution may stain exposed skin or clothing.
2. Eye protection, such as goggles or face shields, should always be used when handling CIDEX OPA.
3. Fluid repellent gown or apron should always be used to prevent OPA from getting on clothing or skin.
4. Use CIDEX OPA Solution in a dedicated room with adequate ventilation and in closed containers with tight-fitting lids. Please reference local/state/country requirements for appropriate number of air exchanges per hour. If adequate ventilation is not provided by the existing air handling system, use in local exhaust hoods, or in ductless fume hoods/portable ventilation devices which contain filter media which absorb *ortho*-phthalaldehyde from the air.

5. Contaminated reusable devices **MUST BE THOROUGHLY CLEANED** prior to disinfection, since residual contamination with soil or lubricants will decrease the effectiveness of the disinfectant.
6. The user **MUST** adhere to the Directions for Use, as modification to the Directions for Use may affect the safety and effectiveness of the disinfectant.
7. Do not use CIDEX OPA Solution on medical devices that are intended for use in a sterile area of the body (e.g. cataract surgical instruments).
8. The reusable device manufacturer should provide the user with a validated reprocessing procedure for that device using CIDEX OPA Solution.
9. The use of CIDEX OPA Solution in automated endoscope reprocessors must be part of a validated reprocessing procedure. Use CIDEX OPA Solution Test Strips to detect *ortho*-Phthalaldehyde concentration before each cycle to detect the MEC. Follow the Directions for Use provided with the CIDEX OPA Solution Test Strips.
10. If using filtering face pieces, please ensure that they are capable of preventing inhalation of vapor. Do not use filtering face pieces that are only designed to prevent inhalation of particulates (e.g. do not use N95 filtering face piece respirator).

DIRECTIONS FOR USE

Cleaning/Decontamination: Blood, other body fluids, and lubricants must be thoroughly cleaned from the surfaces and lumens of medical devices before reprocessing in the disinfectant. Blood and other body fluids should be disposed of according to all applicable regulations for infectious waste disposal.

Refer to the reusable device manufacturer's labeling for instructions on disassembly, decontamination, cleaning and leak testing of their equipment.

Before immersion in CIDEX OPA Solution, thoroughly clean devices, including all lumens, using a cleaning protocol or standard.

Thoroughly rinse and rough dry all surfaces and lumens of cleaned devices.

Usage: **NO ACTIVATION IS REQUIRED.**

Record the date the container was opened on the container label, or in a logbook. After opening, the solution remaining in the container may be stored for up to 75 days (providing the 75 days does not extend past the expiration date on the container) until used. Record the date the solution was poured out of the original container into a secondary container in a log book (separate from the one mentioned above), or on a label affixed to the secondary container. The solution in the secondary container can be used for a period up to 14 days. The product must be discarded after 14 days even if the CIDEX OPA Solution Test Strip indicates a concentration above the MEC.

A. High Level Disinfection: Immerse device completely, filling all lumens and eliminating air pockets, in CIDEX OPA Solution for a minimum of 5 minutes at 20°C (68°F) or higher to destroy all pathogenic microorganisms, including *Mycobacterium bovis*, *Pseudomonas aeruginosa*, pathogenic fungi, and viruses (Poliovirus Type 1; Adenovirus Type 2; Herpes simplex Types 1,2; HIV-1; Influenza Type A [Hong Kong]; Vaccinia; Coronavirus; Coxsackievirus Type B-3; Cytomegalovirus; Rhinovirus Type 42). Remove device from the solution and rinse thoroughly following the rinsing instructions below.

B. Rinsing Instructions

1. RINSING PROCEDURE

a) Manual Processing:

- Following removal from CIDEX OPA Solution, thoroughly rinse the medical device by immersing it completely in a large volume (e.g. 8 litres) of water. Use sterile water unless potable water is acceptable. See item 2 or 3 below.
- Keep the device totally immersed for a minimum of 1 minute in duration, unless a longer time is specified by the reusable device manufacturer.
- Manually flush all lumens with large volumes (not less than 100ml) of rinse water unless otherwise noted by the device manufacturer.
- Remove the device and discard the rinse water. Always use fresh volumes of water for each rinse. Do not reuse the water for rinsing or any other purpose.
- Repeat the procedure TWO (2) additional times, for a total of THREE (3) RINSES, with large volumes of fresh water to remove CIDEX OPA Solution residues. Residues may cause serious side effects, including anaphylaxis in rare instances. SEE WARNINGS. THREE (3) SEPARATE, LARGE VOLUME WATER IMMERSION RINSES ARE REQUIRED.
- Refer to the reusable medical device manufacturer's labelling for additional rinsing instructions.

b) Automated Processing:

- Select a rinse cycle on an automatic endoscope reprocessor that has been validated for use with this product.
- Ensure that the automated rinse cycle selected will thoroughly rinse the medical device including all lumens with large volumes of sterile or potable water equivalent to the reusable device manufacturer's recommendations.
- Verify that each rinse is a minimum of 1 minute in duration unless the reusable device manufacturer specifies a longer time. Ensure that a fresh volume of water is used for each rinse. Do not reuse the water for rinsing or any other purpose.
- Refer to the reusable device manufacturer's labelling for additional rinsing instructions.

2. STERILE WATER RINSE: There are certain devices and procedures that may require sterile water rinsing after high level disinfection. The following devices should be rinsed with sterile water using clean technique:

- Devices intended for use in normally sterile areas of the body.
- Devices intended for use in known immune compromised patients, or potentially immune compromised patients based on institutional procedures (e.g., high risk population served).
- Bronchoscopes, and other semi critical devices may be at risk of contamination when rinsing with potable water. These devices should undergo a sterile water rinse. There may be microorganisms present in the potable water that may be pathogenic for immune compromised patients which may place them at high risk for infection.
- Other semi-critical, lumened instruments should be rinsed with sterile water using a sterile syringe using clean technique.

3. POTABLE WATER RINSE:

- For all other devices, a sterile water rinse is recommended when practical. Otherwise, potable tap water rinse is acceptable.

- When using potable water for rinsing, the user should be aware of the increased risk of recontaminating the device or medical equipment with microorganisms, which may be present in potable water supplies.
- Water treatment systems, such as softeners or deionizers, may add microorganisms to the treated water to the extent that microbial content of the water at the point of use could exceed that of the pretreated drinking water. To ensure proper water quality, adherence to maintenance of the water treatment system(s) is recommended.
- The use of a bacterial retentive (0.2 micron) filter system may eliminate or greatly reduce the amount of these water borne bacteria from the potable water source. Contact the manufacturer of the filter or UV system for instructions on preventative maintenance and periodic replacement of the filter to avoid colonization or formation of biofilms in the filter.
- A device that is not completely dried provides an ideal situation for rapid colonization of bacteria. As these waterborne bacteria are highly resistant to drying, rapid drying will avoid possible colonization but may not result in a device free from these bacteria. A final rinse using a 70% isopropyl alcohol solution can be used to speed the drying process and reduce the numbers of any organism present as a result of rinsing with potable water.

C. Special Instructions for Transesophageal Echocardiography (TEE) probe reprocessing: As with all devices, carefully follow all probe manufacturer recommendations such as use of a sterile protective sheath when performing TEE. Soaking for a minimum of 5 minutes in CIDEX OPA Solution is required for high level disinfection (HLD). Excessive soaking of the probes (e.g., longer than an hour) during HLD and/or not rinsing three times with a fresh quantity of water each time as described in Part B, may result in residual CIDEX OPA Solution remaining on the device, the use of which may cause staining, irritation or chemical burns of the mouth, throat, oesophagus and stomach.

D. Reusage for Disinfection: CIDEX OPA Solution has demonstrated efficacy in the presence of organic soil contamination and microbiological burden during reuse. The *ortho*-phthalaldehyde concentration of CIDEX OPA Solution during its use-life must be verified by the CIDEX OPA Solution Test Strips prior to each use, to determine that the MEC of 0.3% is present. CIDEX OPA Solution may be used and reused within the limitations indicated above for up to 14 days. CIDEX OPA Solution must be discarded after 14 days, even if the CIDEX OPA Solution Test Strip indicates a concentration above the MEC.

MONITORING OF DISINFECTANT:

- During reuse, it is recommended that the CIDEX OPA Solution be tested with the CIDEX OPA Solution Test Strips prior to each use. This is to ensure that the appropriate concentration of *ortho*-phthalaldehyde is present.
- During the usage of CIDEX OPA Solution as a high level disinfectant, it is recommended that a thermometer and timer be utilized to ensure that the optimum conditions are met.
- Visually inspect the solution during the reuse life for the presence of precipitates, which may result from the use of hard water. Discard solution if precipitation occurs.

POST-PROCESSING HANDLING AND STORAGE OF REUSABLE DEVICES:

Disinfected reusable devices are either to be immediately used, or stored in a manner to minimize recontamination. Refer to the reusable device manufacturer's labeling for additional storage and/or handling instructions.

STORAGE CONDITIONS AND EXPIRATION DATE

1. CIDEX OPA Solution should be stored in its original sealed container at controlled room temperature 15 - 30°C (59 - 86°F) in a well-ventilated, low-traffic area.
2. Once opened, the unused portion of the solution may be stored in the original container for up to 75 days until used.
3. The expiration date of the CIDEX OPA Solution is found on the immediate container.

EMERGENCY AND TECHNICAL PRODUCT INFORMATION

For further hazard information please refer to the Material Safety Data Sheet.

Emergency, safety, or technical information about CIDEX OPA Solution can be obtained by contacting your local Advanced Sterilization Products sales representative.

USER TRAINING

The user should be adequately trained in the decontamination and disinfection of medical devices and the handling of liquid chemical disinfectants. Additional information about CIDEX OPA Solution can be obtained by contacting your local Advanced Sterilization Products sales representative.

DISINFECTANT/CONTAINER DISPOSAL INFORMATION

Check local disposal regulations. Glycine (free base) may be used as a neutralizer for CIDEX OPA Solution prior to disposal. A minimum of 25 grams of glycine (free base) should be used to neutralize 3.78 litres of CIDEX OPA Solution. The minimum recommended neutralization time is one hour. Dispose of contents/container in accordance with local/regional/national/international regulations. Do not reuse empty container.

HOW SUPPLIED

Product Code	Description	Case Contains
20391	One Gallon (3.785L)	4 x 3.785L/case Container
20392	CIDEX OPA Solution Test Strips	60 strips/btl; 2 btl/case
20393	CIDEX OPA Solution Test Strips	15 strips/btl; 2 btl/case

For a glossary of symbols used, go to eifu.asp.com.

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© ASP 2021
ASP International GmbH