



GIMA

CLASS 4 MULTI-VARIABLE STEAM STRIPS - dual

Code: 35954
Category: Sterilization Tests
Unit of sale: box of 250 pcs.
Minimum order: 1
Type: No medical device



EAN13: 8023279359541

Description: TYPE 4 MULTI-VARIABLE STEAM STRIPS - box of 250 dual strips

Type 4 multi-variable steam indicator strips are designed to be used in steam sterilizers operating at 134° C - 3.5 min.

The indicator strips will provide assurance of steam penetration into the packs.

The water based, non-toxic chemical indicator in compliance with ISO 11140-1 will show a clear colour change from blue to black after successful steam sterilization.

Test indicator strips are cost effective and can be divided into two separate pieces for economical use.