



GIMA

"ASPEED" SUCTION ASPIRATOR - 230V double pump

Code: 28245
Category: Surgical aspirators and breast pump
Unit of sale: 1 pc.
Minimum order: 1
Type: Medical device
Class: II A
NSIS: 1758504
CND: Z120105
EAN13: 8033717010767



Description: Piston-type continuous cycle electric aspirators give high performance and great durability.
Equipped with a protective thermal cut-out relay.
They require no maintenance or lubrication.
A motor-protector cap totally prevents aspirated liquids or secretions from reaching and damaging the vacuum pump.
Ideal for clinical use.
Made in Italy.

- Aspeed 22 l
- Power: 230 V
- Pump: double
- Case: metal

Technical Specifications: Operating voltage: 230 V-50 Hz
Bottle capacity: 1 l
High vacuum: low flow
Air flow: 22 l/min
Adjustable vacuum level: 0÷ -0.85 bar (0÷ -85 kPa)
Weight: 4.2 kg
Case material: metal
Noise level: 63 dBA

Standard accessories: Autoclavable polycarbonate jar 1,000 cc with safety valve (overflow protection)
Disposable suction liner 1 l
Antibacterial hydrophobic filter
Sterile disposable cannula
Sterile manual flow regulator
Set of atoxic sterilizable silicone tubes
Power Cable
User Manual: GB, FR, IT, DE, ES