



GIMA

SOEHNLE 8352 MULTINA DIGITAL BABY AND CHILD SCALE

Code: 27306
Category: Medical baby scales
Unit of sale: 1 pc.
Minimum order: 1
Type: No medical device



EAN13: 8023279273069

Description: SOEHNLE 8352 MULTINA DIGITAL BABY AND CHILD SCALE

Well-designed, easy to handle baby scale made out of easy to clean plastic.

- suitable to weigh children up to 20 kg removing baby tray
- easy to read LCD-display with 13 mm figures
- simple operation with one button
- fully-electronic scale: no wear and tear parts, no maintenance
- tare function to calibrate for baby blankets/textiles
- hold function: weight value will be held for 10 seconds after unloading
- automatic power-off to save battery capacity

Removable baby weighing tray.

Manual in: DE, GB, FR, IT, ES, NL, DK, SE, HU, FI, PT.

Technical Specifications:

- Capacity: 20 kg
- Dimensions: 601 x 346 x 92 mm
- Battery operated: 9 V
- Sensitivity: 10 g
- Weight: 3.5 kg
- Warranty: 2 years
- Made in Germany