



GIMA

BOWIE DICK TEST PACK

Code: 35950
Category: Sterilization Tests
Unit of sale: box of 20 pcs.
Minimum order: 1
Type: No medical device



EAN13: 8023279359503

Description: BOWIE DICK TEST PACK

- box of 20 pcs
- dimensions: 123 x 123 x 16 mm

Used for daily monitoring of pre-vacuum steam sterilizer.

Based upon the Bowie Dick test result the steam sterilizer is released for usage during the day.

The Bowie Dick test detect failures like air leaks, insufficient vacuum, inadequate steam penetration.

Package and manual in: GB, IT, ES, FR, PT, GR, DE.

- fulfil requirements of EN 285 and ISO 17665
- conform to ISO 11140 - 4 Type 2 / EN 867 - 4 class B (EU standard 7 kilo)
- independently tested to show performance similar to original pack
- clear blue - violet colour change
- easy to interpret and document
- non-Toxic / lead free