



### Product description

Order number	15.1674
Product description	Multivette® 600 Serum Gel CAT, venous and capillary blood samples, preparation: clotting activator/gel, nominal volume: 600 µl, (LxØ) with cap: 53.3 x 8 mm, with print, cap brown, screw cap, flat base, 100 piece(s)/case

### Product characteristics

Type of collection	Venous and capillary
Type of preparation	clotting activator/gel
Colour code	EU/ISO
Label/ Print	with print
Colour of print/label	brown
Cap	standard screw cap
Closure type	screw cap
Application	venous and capillary blood samples
Base shape	flat base
Recommended centrifugation	10000 x g
Centrifugation time	5 Min

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





**Product Specification**  
**Multivette® 600 Serum Gel CAT, 600 µl, cap brown, screw cap**

**Page 2**

**Size**

Sample volume	600 µl
Diameter	8 mm
Length including cap	53.3 mm
Length excluding cap	49.8 mm

**Material & colours**

Product material	Polypropylene (PP)
Colour of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	brown

**Purity & certification**

Product category	In-vitro diagnostic, CE
CE certified	CE
Batched	yes

**Packaging**

Minimum order qty.	1000
Type of smallest subpackaging	case
Piece(s) / inner box	100
Piece (s) / outer case	1000
Piece(s) / pallet	63000
Depth of box	201 mm
Width of box	81 mm
Height of box	81 mm
Depth of case	423 mm
Width of case	213 mm
Height of case	185 mm
Case volume	0.0167 cbm
Weight of product	0.0024 kg
Weight of case	2.97 kg
EAN of inner box	4038917087421
EAN case	4038917005852

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:  
2024-02-21

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**