

# evercare<sup>®</sup> Medical face mask with earloops, type II

REF: 2882

**evercare<sup>®</sup> Medical mask with earloops, Type II** provides a barrier to minimize the direct transmission of infective agents between staff and patient.

## Key Features:

- ▶ A type II medical mask with earloops
- ▶ Bacterial Filtration Efficiency (BFE) ≥ 99%
- ▶ Anatomic pleated folding limiting inward leakage
- ▶ Flexible nose wire for proper fit
- ▶ Color: White



REF	DESCRIPTION	SIZE mm	PACK COUNT dispenser/carton
2882	evercare <sup>®</sup> Medical face mask with earloops, white	165 x 90	50/2500

**Sales Unit:** Transport Carton

# Technical Summary

## Regulations

The product complies with legal requirements for medical devices



**REACH regulation:** The product(s) covered by this data sheet does(do) not contain substances classified as carcinogenic, mutagenic or toxic for reproduction (CMR), including any of the Substances of Very High Concern (SVHC) as listed in the latest available version of the Candidate List published by European Chemicals Agency (ECHA).

**Biological evaluation:** EN ISO 10993-5, EN ISO 10993-10

**Compliance with legal requirements:** Medical Devices Directive 93/42/EEC, Class I

**Compliance with product standards:** EN 14683 Type II, EN ISO 11737-1

**Quality standards:** EN ISO 13485

**Label information and packaging compliance:** EN ISO 15223-1, EN 1041

## Material

Product Part	Raw Material & Characteristic
Outer coverstock	Spunbond PP nonwoven, white
Filter media	Meltblown PP nonwoven, white
Inner coverstock	Spunbond PP nonwoven, white
Nose wire	Plastic
Ear loops	Spandex, Polyester and Polyamide, white, latex-free

## Warnings and Recommendations

**Warnings:** The product is single use. If reused performance of product may deteriorate, cross contamination may occur.

**Storage Recommendation** Recommended storage in clean, dry space in ambient temperature. Products should be protected from direct sunlight, other intensive light sources and ozone.

**Disposal Recommendation** Product and packaging materials to energy recovery or mixed waste.

**Shelf Life** 3 Years

## Packing Information

Package	Material	Size ( L x W x H ) mm
Single Pack	N/A	N/A
Dispenser Pack	Inner PE bag, LDPE; Carton, TCF quality	188 x 98 x 87
Transport Carton	Carton, TCF quality	460 x 395 x 515
Packing Tape	Transparent, PVC-free	N/A

## Barcodes

REF	Single Pack	Dispenser	Transport Carton
2882	N/A	 6 438129 324624	 6 438129 524628



Manufactured in China for  
**OneMed Group Oy**  
 Address: Metsäläntie 20  
 Postal Address: P.O. Box 10  
 FI-00321 Helsinki, Finland  
 Phone +358 20 786 6800  
 www.onemed.com  
 Domicile Helsinki  
 Business ID 2039640-1