

# THERMOMETERS

DIGIVET, DIGIFLASH, DIGIFIX



injection



surgery



orthopedics



hospitalization



capture



protection



hygiene



livestock



grooming

**génia**  
the veterinary preference

## DIGIFLASH

Digital thermometer for pets and large animals.

- **Flexible** tip, length 7 cm.
- Temperature measurement in 5 seconds.
  - Recording the last 5 values.
  - Temperature range: 32 to 43 °C / 90 to 110 °F.
  - Accuracy: 0.1 °C / 0.18 °F.
  - Alarm indicating the end of measurement.
- Low battery symbol display.
- Auto power-off.
- Waterproof.
- Shock resistant.

Supplied in a transparent polycarbonate box with an additional battery.



FAST SETTING | FLEXIBLE  
5 SECONDS | TIP

Code	Length	Tip length
227972	15,5 cm	7 cm

**Outer packaging :**  
Master case of 12 and 240 units.

## DIGIFIX

Digital thermometer for pets and large animals.

- **Rigid** tip, length 7,5 cm.
- Temperature measurement in 10 seconds.
  - Recording the last value.
  - Temperature range: 32 to 43°C / 90 to 110°F.
  - Accuracy: 0,1°C / 0,18°F.
  - Alarm indicating the end of measurement.
- Low battery symbol display.
- Auto power-off.
- Waterproof.
- Shock resistant.

Supplied in a transparent polycarbonate box with an additional battery.



RIGID | TIP

Code	Length	Tip length
227978	15,5 cm	7,5 cm

**Outer packaging :**  
Master case of 12 and 240 units.

## DIGIVET

Digital thermometer for pets.

- **Flexible** tip, length 5 cm.
- Temperature measurement in less than 10 seconds.
  - Temperature range: 32 to 42 °C / 90 to 108 °F.
  - Accuracy: 0.1 °C / 0.18 °F.
  - Alarm indicating the end of measurement.
- Low battery symbol display.
- Auto power-off.
- Waterproof.
- Shock resistant.

Supplied in a transparent polycarbonate box with an additional battery.



FLEXIBLE | FAST SETTING  
TIP | 10 SECONDS

Code	Length	Tip length
227970	12,5 cm	5 cm

**Outer packaging :**  
Master case of 12 and 240 units.



HOSPITALIZATION