



## Welch Allyn MacroView™ Veterinary Otoscope

FOR VETERINARIANS



### A Difference You've Got to See to Believe

Get a wider, deeper, clearer view of the ear with the Welch Allyn MacroView™ Veterinary Otoscope, the latest advance to the device that helps you perform one of your most common procedures. Using advanced technology in an innovative design, the MacroView™ Veterinary Otoscope gives you a large, sharp, nearly full view of the tympanic membrane—even in large dogs—virtually eliminating the need for panning. And with an easy-to-focus eyepiece, you can even adjust focus for different animals and ear canal lengths, meaning you get crystal-clear views of the ear canal at any point—allowing you to more effectively diagnose and treat your patients faster than ever before.

### The Welch Allyn MacroView™ Veterinary Otoscope gives you:

- More detail with better clarity—approximately 27% greater magnification than with traditional otoscopes
- Approximately 35% greater field of view than with traditional otoscopes
- The ability to adjust focus for variable ear canal lengths and can correct for farsighted eyes
- Tip Grip for secure fastening of ear specula, along with ejection feature for easy removal
- Built-in throat illuminator


**WelchAllyn®**

Advancing Frontline Care™

# The Welch Allyn MacroView™ Veterinary Otoscope Gives You a Bigger, Better, Sharper View



## INCREASED FIELD OF VIEW




**MacroView™  
Veterinary Otoscope**

Get a nearly full, crystal-clear view of the tympanic membrane—even in large canines!


*"The wider angle of view and ability to focus the MacroView™ Veterinary Otoscope means better evaluation by my doctors and better outcomes for our patients. Finally, I'm able to fully visualize the ear canal and tympanum of my pet patients!"*

—Dr. Ernest Ward, Jr.  
Seaside Animal Care, Calabash, N.C.



**Traditional Otoscope**

## LONGER, ADJUSTABLE FOCAL LENGTH



## Ordering Information

23862 MacroView Veterinary Otoscope with Specula

### Accessories

23834 4 mm MacroView Veterinary Reusable Specula  
 23835 5 mm MacroView Veterinary Reusable Specula  
 23837 7 mm MacroView Veterinary Reusable Specula  
 23842 MacroView Veterinary Reusable Specula Set  
 23847 MacroView Veterinary 7 mm Instrumentation Specula  
 06500-U Replacement Lamp

*Diagnostic Sets featuring the MacroView Veterinary Otoscope are also available.*

To find out how the new Welch Allyn MacroView™ Veterinary Otoscope can give you bigger, sharper views of the ear, go to [www.welchallyn.com/dvm](http://www.welchallyn.com/dvm) or call 1.800.535.6663.

4341 State Street Road, Skaneateles Falls, NY 13153 USA  
 (p) 1.800.535.6663 (f) 1.315.685.3361 [www.welchallyn.com/dvm](http://www.welchallyn.com/dvm)

**WelchAllyn®**

Advancing Frontline Care™