

# SPLINTS RANGE

FLAT, ROLL, CONCAVE SPLINTS



injection




surgery



orthopedics




hospitalization



capture



protection



hygiene



livestock



grooming

**génia**  
the veterinarian's preference

## SPLINTS

Concave plastic harnesses; red, blue or green depending on size; coated with a thick protective foam.

- Very strong material, poly-gloss finish.
- No sharp angles, thick closed-cell protective foam that extends beyond the edges of the splint to protect the animal.
- Can be cut to adjust to the animal's member size.
- Air holes along the length of the splint allow air to circulate and prevent scab formation.

Supplied by unit or 12.



NO SHARP ANGLE | VERY STURDY MATERIAL

Code x 1	Code x 12	Length	O/pack.
307801	307800	20 cm	x 24
307806	307805	30 cm	x 8
307811	307810	40 cm	x 6

## FLAT SPLINTS

An aluminium side and a protective soft foam side.

- Easy to cut with Mayo serrated scissors to adjust to the animal's limb.
- Radiotransparent.
- Sturdy, light and atraumatic.
- Ideal for a temporary immobilization.
- Available in 5 widths.

Supplied by 12.



RADIOTRANSSPARENT | ATRAUMATIC

Code x 12	Size
307820	2 x 20 cm
307822	2,5 x 30 cm
307824	5 x 40 cm
307826	8 x 45 cm
307828	10 x 45 cm

## ROLL SPLINTS

NEW

Aluminium roll splint with foam on both sides.


- Easy to cut.
- Cut with Mayo serrated scissors to fit the animal's limb.
- Radiolucent.
- Flexible, atraumatic, easily takes the shape of the limb requiring immobilisation.
- Available in 3 widths.

Supplied by 4.



SOFT | RADIOTRANSSPARENT | ATRAUMATIC

Code x 4	Size
307830	4,5 x 25 cm
307832	5,5 x 35 cm
307834	9 x 40 cm



HOSPITALIZATION