



# Midmark B60 SI Sterilizer



## SAVE TIME WITH FASTER CYCLES

Larger capacity and faster cycle times help you maximize the efficiency in your practice. With a built-in Steam Generator, the Midmark DentiSure B60 SI completes most cycles in just one hour, and the most demanding cycles in 90 minutes.

## ENJOY HIGH-PERFORMANCE AND RELIABILITY

Save money by enjoying lower water consumption levels, using only 1.5L per cycle, easily measure water quality, and see perfect-dry results. Reverse Osmosis cartridges now have a longer life, so your expenses are significantly reduced.

## EASY TO USE

Our most advanced user interface is quick to understand and very easy to use. Features allow you to schedule test cycles with “start delay” and save energy with a quick setting that stops energy consumption when the autoclave is not in use. The unit is not only easy to use, but easy to service when necessary with built-in RO cartridges and easily removable panels.

- 
- Available with or without cabinet.



## SPECIFICATIONS

SPECIFICATIONS	B60 SI
Chamber material	Stainless steel
Overall Width (mm)	680
Overall Height (mm) (with desktop)	820 (1600)
Overall Depth (mm)	783
Weight (kg) (with desktop)	146 (204)
Number of Trays	4
Tray Width (mm)	310 / 220
Tray Depth (mm)	550
Tray Height (mm)	25
Supply Voltage (VAC +/- 10%)	400 VAC (3 phases)
Mains Frequency (Hz)	50-60
Class	B
Solid Loads (unwrapped)	Yes
Wrapped Materials	Yes
Narrow Lumen	Yes
Simple Hollow Items	Yes
Porous Loads	Yes
Small Porous Loads	Yes
Number of Preprogrammed Cycles	7
Helix and Vacuum Test	Yes
Bowie&Dick e Vacuum Test	Yes
Separate Clean and Dirty Water Tanks	Yes
Built in Water Conductivity Meter (standard)	Yes
Instant Vaporizer/Steam Generator	Yes
Door Closure (motor-operated)	Double-motors
Printer (standard)	Yes
USB integrated (optional)	Yes
USB external (optional)	Yes
Graphic Display	Yes
Night Cycle	Yes
MTS (Last 10 Cycles Reprint)	Yes
Parts Warranty	2 years/ 1500 cycles
Chamber Warranty	5 years/ 5000 cycles

STANDARDS: EU Directive 93/42/EEC (medical devices); EU Directive 2006/42/EEC (machines); EU Directive 2014/68/UE (pressure equipment); EU Directive 2014/30/UE (EMC); EU Directive 2014/35/UE (LVD); UNI EN 13060:2015 (Small Steam Sterilizers); CEI EN 61010-1:2013 and CEI EN 61010-2-40:2016 (laboratory equipment)

## Related Products



Also available in our smaller B Class units (B23, B18 and B6 models)



QuickSeal™ Heat Sealer



QuickClean™ Ultrasonic Cleaner

www.newmedsrl.com  
 Newmed srl - a Midmark Company, Via Lenin 79/A  
 42020 Quattro Castella, Reggio Emilia, Italy.  
 Tel: +39 0522 875 166 Fax: +39 0522 243 096  
 newmed.sales@midmark.com  
 ENMI-402 Rev.2 09/19



  
 midmark™