

Surgical Drapes

Deliver industry-leading, all-around performance



**Draped in
confidence**


WHY MEDLINE SURGICAL DRAPES AND GOWNS?

Operating theatres demand best product performance at all times.

That is why Medline focusses on consistently delivering surgical drapes and gowns that are reliable and driven by performance. Listening to the feedback from our customers, we strive to develop products that meet and exceed the needs of operating theatre professionals.

The benefits of Medline surgical drapes and gowns

- » No. 1 supplier of drapes and gowns in the world
- » Used in nearly 23 million surgical procedures across Europe
- » The broadest product line in the industry
- » Developed for clinicians and by clinicians
- » Adapted to customer needs
- » Best-in-class drape and gown fabrics
- » Adaptability of a family-run business but with the scale and structure that comes from being a major healthcare company
- » Maximises value thanks to our distribution capabilities and custom packs



2. PHYSICAL STRENGTH

Heavier-weight and tear-resistant fabrics provide better protection, especially in high-movement cases like orthopaedics.




1. FLUID CONTROL

Highly absorbent fabrics and fluid-impervious zones create a protective barrier and help prevent blood strikethrough and fluid contamination.



3. DRAPEABILITY AND USEABILITY

Innovative fabric technology and design features help ensure patient safety, improve drapeability and save preparation time in the operating theatre.



4. IDEAL FIT

Gowns in soft, breathable fabrics fit comfortably and give wearers complete range of motion.

5. LINT AND ABRASION RESISTANCE

Industry-standard lint- and abrasion-resistant drapes stand up to the rigours of surgery and help reduce the risk of cross-contamination.



INNOVATION ALL AROUND


The operating theatre demands our best performance.

That's why Medline focusses on elevating performance without compromising value. We partner with our customers to design procedure-specific drapes that meet your expectations, then take it a step further.


Designed to excel in all the right places

- » Generous dimensions provide a complete sterile field.
- » Strong hook-and-loop line holders secure lines well and save preparation time.
- » Hand-placed adhesive strip provides exceptional precision and protection at the incision area.
- » Anti-slip instrument pad keeps instruments at hand whilst preventing them from falling to the floor.
- » Critical-zone absorbency pad features the highest absorbency rating and abrasion resistance where it's needed most.*

Nearly 23 million surgical procedures were performed in 2018 across Europe with our products.



*Data on file




OPST[™] ESSENTIAL SURGICAL DRAPES

MADE OF A BILAMINATE, FLUID-ABSORBENT AND IMPERVIOUS MATERIAL

Medline's Essential two-layer material not only absorbs liquids and bodily fluids, but its inner layer also provides an impervious and complete barrier against fluids and bacteria.*

- » Advanced composite technology provides excellent drapeability and patient comfort.
- » Low-linting helps reduce the risk factors associated with infection.
- » Essential is available with and without enhanced high-absorption reinforcement in critical zones.
- » This drape meets or exceeds a high performance level according to EN 13795.*

*Data on file




OPST[™] ULTIMATE SURGICAL DRAPES

A MEDLINE EXCLUSIVE THAT IS FLUID-ABSORBENT AND IMPERVIOUS

Made of a three-layer compound of different materials, the Ultimate drape lives up to its name—absorbing fluids throughout the entire drape, providing 100 per cent barrier strength against fluids and bacteria and ensuring patient comfort.

- » Ultra-low linting helps reduce the risk factors associated with infection.
- » Fluid-impervious and absorbent trilaminar material responds to the surgical team's needs in the operating theatre.*
- » A soft non-woven layer helps manage a patient's moisture during long surgeries.
- » This drape has a high performance level according to EN 13795.*

*Data on file




OPS™ ADVANCED SURGICAL DRAPES

THE INDUSTRY STANDARD WITH 5-LAYER STRENGTH

Composed of five strong fibre layers, our SMMMS drape provides drapeability, softness and a barrier against micro-organisms.*

- » Latest SMS fabric technology ensures exceptionally high barrier strength.
- » The 5-layer microporous structure offers excellent breathability and helps manage a patient's moisture.
- » Low-linting helps reduce the risk factors associated with infection.
- » Heavy-duty absorbent reinforcement in critical zones helps prevent blood strikethrough and fluid contamination.

*Data on file




ECLIPSE SURGICAL DRAPES

FABRIC WITH A COTTON-LIKE APPEARANCE AND FEEL

Medline's surgical drapes' fabrics meet or exceed the industry standards.

- » Eclipse is drapeable, repellent and breathable, and offers textile-like qualities.
- » Low-linting helps reduce the risk factors associated with infection.
- » Eclipse can be used with or without an additional absorbent reinforcement in critical zones.
- » This drape meets or exceeds the industry drape performance requirements.



EXCLUSIVE FEATURES OF OUR DRAPES AND PACKS PORTFOLIO

1. PREMIUM MEDICAL ADHESIVE TAPE



2. TRACEABILITY STICKERS

4 traceability stickers
available on every product
pouch



3. CLEAR LABELS

Clear product labels with item
number, description and large
drawings



4. COLOUR CODING

Per pack for differentiation
between materials



5. REDUCE WASTE

Medline man to show you the
patient's position and how
to unfold the drape, thereby
reducing improper unfolding
and drapes that have to be
discarded



6. SMS WRAP

On all our products



**Choose your
protection**



**ALWAYS
ON.**

Medline Industries Ltd
3rd Floor
Quayside Wilderspool Business Park
Greenalls Avenue
Warrington WA4 6HL
United Kingdom
Tel.: +44 844 334 5237
Fax: +44 844 334 5238

www.medline.eu/uk
uk-customerservice@medline.com

These products are class I medical devices intended to be used by healthcare professionals. Before use, consult instructions and precautions on the corresponding labelling.

We reserve the right to correct any errors that may occur within this brochure.
© 2020 Medline Industries, Inc. Medline is a registered trademark of Medline Industries, Inc. ML610-EN 01/2020.



Medline International France SAS
5 Rue Charles Lindbergh
44110 Châteaubriant, France

CE 2797 (BSI - The Netherlands)