

# Veterinary Fluid Warmer



# Veterinary Fluid Warmer



## FW2-VET

Fluid warmer is a device heating fluid inside infusion set /blood transfusion set based on thermal transfer principle. It accepts various brand of infusion set/blood transfusion set.

### Fluid Warmer Features

- Compact and light weight
- Colorful LCD screen with detailed info
- Double groove heating for two infusion sets
- Optional working with drop sensor
- Fast heating, 3~5 minutes to reach the set temperature
- With over-temperature & low temperature alarm
- Cut off the power supply automatically when overheat of fluid
- Two flexible installations: with hanger or pole clamp
- The pole clamp has hidden function, saving space

### Specification

Tube size	3~4.5 mm
Heating range	30 °C ~ 50 °C, with increment of 1°C or 0.1 °C
Accuracy	≤ ± 1 °C
Classification	Class I / type B
Power	AC 100~240V, 50Hz/60Hz
Net weight	550g
Packing dimension:	22.0*19.0*12.1 cm

