

# PrO-VetPCVs



Display indicative only.

**Dims HWD** .....235 x 235 x 350mm

**Weight**.....12.8 Kg (without rotor)

**Power**.....125 watts

**Rotor**.....BRK5401 & BRK5408

Program	Speed	Time
<b>Program 5401</b>		
Haematocrit (Pcv)	12,000Rpm	300Secs
Serum 1	6,000Rpm	600Secs
Serum 2	5,000Rpm	300Secs
Urine	3,500Rpm	300Secs
<b>Program 5408</b>		
Serum 1	6,000Rpm	180Secs
Serum 2	6,000Rpm	240Secs
Serum 3	5,000Rpm	300Secs
Urine	3,500Rpm	300Secs

## Serum (15ml max) & Haematocrit Rotor



### Rotors

	BRK5401	BRK5408
Rotor type	24 x Capillary & 12 x 2ml	8 x 15ml
Tube size max	2 x 75mm & 11 x 40mm	17 x 120mm
Minimum speed	500 Rpm	500 Rpm
Maximum speed	12,000 Rpm	6000 Rpm
Maximum Rcf (G)	13500	3500
Radius max cms	8.5	8.50
Sample tube angle (α)	0 & 60	30
Acceleration time (secs)	30	30
Deceleration time (secs)	30	30
Autoclaveable (frequency)	121°C (10)	121°C (20)

### Accessories



**8003 PCV reader**



**8014 Sealing clay**

Part Number	Description
8012	Plain Capillary tubes pack 1000
8013	Heperanised Capillary tubes pack 1000