

Next2Vet

HydroJet 2

Ny uppgraderad hydrojet från HDE. Effektiv mun skölj mellan tänder och i foderfickor



Hydro jet II Artikel no 709050

The Hydro jet II includes the following:

Hydro Jet II with battery and charger

Rinsing nozzle

Diastema nozzle

6-meter water tube

1-liter bottle

Slim carriage bag

HDE NO: ACC_SHJ2

The Hydro jet II has been completely redesigned for easier work with several new smart features:



Fast & easy brass nozzle connectors



Water flow regulation from 0,75 to 1,5 liter/minute



Battery level indicator



Curved rinsing nozzle



Rounded-off diastema tip



Battery moved on the top

Waterpressure: 7 bars

Battery usage time: 2 hours of continuous use

Battery charging time: 1,5 - 2 hours

The diastema nozzle is perfect for flushing diastemas and removing interdental food residues. Its length and shape allow easy access in the back of the mouth.

The rinsing nozzle is for general mouth cleaning all along the different dental treatment stages. Its specific curved shape naturally fits the mouth.

För mer info kontakta oss på Next2Vet

Next2Vet

info@next2vet.se

0411-55 57 30