



# VetPro<sup>®</sup> 1000

Mobile Dental  
Delivery System



# The VetPro® 1000 Mobile Dental Delivery System

## Oil-Free, Worry-Free, Hassle-Free

A reliable veterinary dental delivery solution that helps you provide better patient care.



### Oil-free compressor

Our exclusive oil-free compressor technology allows you to provide your patients with better, more efficient care without the mess or hassles of oil-lubricated units. The VetPro® 1000 dental delivery system requires no oil changes and very little maintenance for hassle-free operation. Oil-free means there is no oil to contaminate your handpieces, providing a clean, dry, aseptic work environment for you, your staff and your patients.

### Designed for veterinary dentistry

The VetPro® 1000 has precise controls designed specifically for veterinary dentistry. You can switch from handpiece to handpiece without touching a button – the unit will automatically adjust as the handpiece is lifted from its cradle. The oil-free compressor means there are no oily parts for animal hair to stick to, and, it's one of the quietest compressors on the market with a sound level of 40 dB.

### Durable construction

Its durable construction withstands the rigors often found in clinical settings. The VetPro® 1000 is a completely enclosed system to protect internal mechanics from contaminants and animal hair, and is made from sturdy ABS plastic to withstand the strains caused by frequent usage.

### Supported by Midmark

The VetPro® 1000 gives you unsurpassed reliability, backed by our trained service technicians and independent service network. You will have peace of mind knowing that over 90 years of manufacturing experience goes into each product. VetPro® 1000...the oil-free, worry-free, hassle-free solution for your veterinary dental needs.

#### Specifications

Overall depth: 18.5" (46.99 cm)  
Overall width: 23" (58.42 cm)  
Overall height: 36.5" (92.71 cm)  
Casters: 3" (7.62 cm)  
medical grade

#### Compressor

Oil-free 1/3 hp motor  
Working pressure: 70-100 psi  
Flow rate: 2.5 CFM at 40 psi  
Sound level: 40 dB  
Air reservoir: 1 gallon  
Manual air tank purge

#### Water Module

Two 1-liter reservoirs (toggled)  
Pressurized system: 30-40 psi  
Aseptic, non-retracting  
HPC coolant  
Manual handpiece flush system  
External quick connects for  
water and air

#### Electrical

Voltage: 115 VAC, 6 Amps  
Frequency: 50/60 Hz  
UL 60601-1:2003  
CAN/CSA-C22.2 No. 601-1-M90  
Powercord: 10' hospital grade

#### Accessories

VetPro® Lube-Free High-Speed and Low-Speed Handpieces included  
Optional, Integrated Acteon Satelec® SP Newtron® Piezo Scaler, with or without LED Light  
Handpiece connection type: Midwest 4-hole style  
Optional fiber optics  
Bur starter kit included  
10 Disposable prophylaxis angles included

#### Model Numbers

8000-001	VetPro® 1000 with integrated scaler
8000-002	VetPro® 1000 without scaler
8000-003	VetPro® 1000 with integrated scaler and fiber optics
8000-004	VetPro® 1000 without integrated scaler; with fiber optics
8000-007	VetPro® 1000 with LED scaler and dual fiber optic high-speed handpieces
8000-008	VetPro® 1000 with LED scaler and fiber optics



facebook.com/MidmarkAnimalHealth

midmarkanimalhealth.com

Midmark is an ISO 13485 and ISO 9001 Certified Company. • 1-800-MIDMARK • 1-937-526-3662 • midmark.com  
© 2011 Midmark Corporation • Products subject to improvement changes without notice • Litho in U.S.A. • 007-0725-00 Rev. E1 (4/11)



Because we care.