

Veta 5/ Veta 3

Veterinary Anesthesia Machine

Redefining Animal Anesthesia Machine



www.mindrayanimal.com

P/N:ENG-Veta 5/ Veta 3-210285X6P-20210420
©2021 Shenzhen Mindray Animal Medical Technology Co., Ltd. All rights reserved.




Anesthesia safety has always been a key focus for veterinarians and pet owners, but as veterinary medicine advances, we meet new challenges and higher expectations.

The Veta series of anesthesia machines combines advanced clinical practices with cutting-edge technology. Its Intelligent Ventilation mode and integrated design transforms the veterinary anesthetist experience, ensuring the safety of patient and anesthetist alike.

Safety Renovating

Intelligent ventilation mode tailored to the patient



Better display of patient status


- Waveforms and parameters display
- Optional sidestream CO₂ monitoring



More safety features

Auto system leak test

- Illustrated on-screen instructions
- One-step startup



Waste gas canister replacement reminder

- Accurately calculates weight changes
- Automatic canister change reminder

Multiple types of waste gas scavenging system



Gas filter canister



Passive scavenging



Active scavenging

Form Refreshing


Integrated design creates a simple and comfortable operating space



Integrated ventilator-no hoses means lower risk of leaks



8-inch capacitive touch screen, intuitive operation and smoother interaction



Electronically driven by a high-performance turbine, no drive gas needed



Flexible workstation

Multiple mounting solutions for patient monitors, pumps and oxygen generators improves utilization of space.



Modular upgrade

Upgrades from Veta 3 to Veta 5 are available, adding ventilators and providing more options.



Veta 3

Veta 5