

M3B

Vital Signs Monitor



Outpatient



Transport



Ward



Portable
Design



Respiration CO₂



EDAN

A world of potential

M3B Vital Signs Monitor

With SpO₂ (Nellcor Oximax™ optional) and Respironics CO₂ monitoring technologies, the M3B ensures effective capnography monitoring for mechanically ventilated and non-intubated patients during continuous long-term monitoring.

Measured Parameters: SpO₂ + CO₂

- 5.6" high resolution color TFT-LCD screen
- Lightweight & portable design
- User-friendly interface
- Flexible configurations to meet different clinical needs
- Real-time measurements display
- Trend table display
- Built-in rechargeable Lithium-ion battery for 14 hours of continuous work
- USB flash disk support

72 h
Trend Review


800
Alarm Review

- Barcode scanner support
- Bi-directional with MFM-CMS central monitoring system
- Compatible with CMS-Lite for PC data management
- Wired LAN/Wi-Fi support
- HL7 support via XML files or MFM-CMS
- Nurse call
- Suitable for adult, pediatric and neonatal patients

Seamless Connectivity



Mounting Solutions



Specification

Configuration	SpO ₂	Respironics CO ₂ (Sidestream/Mainstream)
M3B(SpO ₂ + CO ₂)	✓	✓
M3B(CO ₂ only)	—	✓



Global Headquarters:

Edan Instruments, Inc. | 15 Jinhui Road, Pingshan District, Shenzhen
518122 P.R. China | +86.755.26898326 | www.edan.com | info@edan.com

U.S. and Canada inquiries:

EDAN Diagnostics, Inc. | 9918 Via Pasar, San Diego, CA 92126
+1.858.750.3066 | www.edandiagnostics.com | edan-info@edandiagnostics.com

© Edan Instruments, Inc. All rights reserved. Features and specifications are subject to change without prior notice. No reproduction, copy or transmission may be made without written permission. Not all products or features are available in all countries, contact Edan for local availability.



ENG-PM-M3B-V2.6-20230516