



GIMA

CHLOREXIDINA SCRUB KIT - Multilanguage


Code: 36574

Category: Brushes

Unit of sale: box of 30 pcs.

Minimum order: 1

Type: PMC - disinfectant



EAN13: 8023279365740

Description: CHLOREXIDINA SCRUB KIT dispenser of 30 pcs - Multilingual
Label: GB, FR, ES, DE, PT, PL, NL, GR (blister and dispenser), PL, RU, Arabic (blister only)

Italian Health Ministry registration PMC 19531

With solution of Chlorhexidine Gluconate 20% 4 g in 100 ml

Sterile hand-cleaning hospital kit.

It includes a scrub, a sponge and a plastic nail cleaner.

Dry, without detergent.

Sterilization method: ETO.

All scrub kits are individually blistered in 30 pieces dispenser

Made in Italy.

All scrub kits are individually blistered in 30 pieces dispenser.