



# GIMA

## BIOPSY PUNCHES diameter 2 mm

**Code:** 22612

**Category:** KAI, GIMA and SIEFEL biopsy punches and dermal curettes

**Unit of sale:** box of 20 pcs.

**Minimum order:** 1


**Type:** Medical device

**Class:** II A

**NSIS:** 1548778

**CND:** A010299

**EAN13:** 4560146922605



**Description:**

- Biopsy Punches Ø 2 mm
- Blade length: 7 mm
- Stainless steel sharp cutting edge, ribbed handle for sure grip.
- Size embossed on each punch for easy identification.
- Suggested for dermatology, gynaecology, podiatry.
- Each punch is individually sealed in a peel-apart see-through package.