





-  **BEST choice for dogs' and cats' care!**
*The one and only one product in the world
with full medical grade supports professional power.*

Easy to select cuffs



Features

- ↘ Three-color LEDs provide better recognition
- ↘ Big Memory size up to 999 sets
- ↘ Cyclic (auto) measurement power up for surgery purpose
- ↘ Large LCD gives easiest user-interface to do setup
- ↘ Rechargeable lithium battery does 800 times of operation
- ↘ The only medical grade NIBP you can hold so easily
- ↘ Full range cuffs for selection



Specifications

Non-Invasive Blood Pressure

Technique	Oscillometric
Parameters	Systolic Pressure, Diastolic Pressure, Mean Arterial Pressure, Heart Rate
Reading Time & Range	<40 seconds; 30 to 300 mmHg
Measurement Cycle	1, 2, 3, 5, 10, 15, 30, 60, 90 minutes, 2 hrs; or off
Resolution & Accuracy	Resolution: 1mmHg; Accuracy: ± 3 mmHg or $\pm 2\%$ of reading, whichever is greater, over full range
STAT mode	5 minutes of continuous reading

Heart (Pulse) Rate

Range	Max: 32 to 240 bpm
Heart Rate Resol & Accuracy	Resolution: 1bpm; Accuracy : ± 2 bpm or $\pm 2\%$

Record

Data Storage	999 measurements
Data output	Mini USB port

Power

Supply Voltage (Recharge)	100 to 240 VAC; 47/63 Hz ; Recharge time: 7hrs
Battery Specification	Lithium battery; Test : up to 600 measurements (1 min Cycle)
Power Save Mode	Automatic power off if idle for 30 minutes

Environment

Operating & Storage Temp	Operating : 5°C to 45°C; 41°F to 113°F Storage : -10°C to 50°C; 14°F to 122°F
Operating Humidity	15% to 90%, non-condensing

Display

Type	LCD / LED
Display Status	Indicate "low battery", "external power"
Alarm Indication	Audible; Visual

Mechanical

Weight & Dimension	475g with battery; 190 mm(L) x 90 mm(W) x 50 mm(H)
--------------------	--

Classification

Type of Protection	Class I equipment
--------------------	-------------------