

# PrO-VetSerum



Display indicative only.

**Dims HWD** .....235 x 235 x 350mm  
**Weight**.....12.8 Kg (without rotor)  
**Power**.....125 watts  
**Rotor**.....BRK5408 & BRK5424

Programs	Speed	Time
Serum 1	6,000Rpm	180Secs
Serum 2	6,000Rpm	240Secs
Serum 3	5,000Rpm	300Secs
Urine	3,500Rpm	300Secs

## Serum Rotor to 15ml *Reducers available for 5 & 10ml*



Rotors	<b>BRK5408</b>
Rotor type	8 x 15ml
Tube size max	17 x 120mm
Minimum speed	500 Rpm
Maximum speed	6000 Rpm
Maximum Rcf (G)	3500
Radius max cms	8.50
Sample tube angle (α)	30
Acceleration time (secs)	30
Deceleration time (secs)	30
Autoclaveable (frequency)	121°C (20)

## Reducers



Rotor	BRK5408
Part number	<b>RM05 (5ml)</b>
Tube size max	13 x 80mm
Part number	<b>RM10 (10ml)</b>
Tube size max	16 x 100mm