Transducers for **Vetus E7 Ultra**



mindray
animal medical

Probe	Illustration	Type of Array	Parameters	Application	Biopsy kit
C5-1s		Curved	• Bandwidth: 1.2 – 6.0 MHz • Extended FOV: 72° • Convex Radius: 60 mm	Abdomen (Canine, Equine, Bovine, Ovine), Reproduction (Canine, Equine, Bovine, Ovine), MSK (Canine, Equine, Bovine, Ovine)	NGB-022, available, multi angle, reusable; LPUBKG60, disposable; CIVCO 658-004, disposable
C6-2Gs		Curved	• Bandwidth: 2.6 – 8.2 MHz • Extended FOV: 104° • Convex Radius: 20 mm	Abdomen (Canine, Equine, Bovine, Ovine), Cardiology (Canine), Reproduction (Canine, Equine, Bovine, Ovine)	NGB-024, available, multi angle, reusable; LPUBKG81, disposable
C11-3s		Curved	• Bandwidth: 2.6 – 12.8 MHz • Extended FOV: 114° • Convex Radius: 15 mm	Abdomen (Canine, Feline), Cardiology (Canine, Feline), Reproduction (Canine, Feline), Small Parts (Canine, Feline), MSK (Canine, Feline)	NGB-018, available, multi angle, reusable
mSC6-1s		Curved	• Bandwidth: 2.6 – 8.2 MHz • Extended FOV: 119° • Convex Radius: 20 mm	Abdomen (Canine, Feline), Cardiology (Canine, Feline),Reproduction (Canine, Feline), Small Parts (Canine, Feline), MSK (Canine, Feline)	NGB-063, available, multi angle, reusable

Probe	Illustration	Type of Array	Parameters	Application	Biopsy kit
C13-3s		Curved	• Bandwidth: 2.6 – 12.8 MHz • Extended FOV: 144° • Convex Radius: 11.5 mm	Abdomen (Canine, Feline), Cardiology (Canine, Feline), Reproduction (Canine, Feline), Small Parts (Canine, Feline), MSK (Canine, Feline)	not available
L13-3Ns		Linear	• Bandwidth: 3.0 – 13.0 MHz • Field of View (max): 3.80 cm	Abdomen (Canine, Feline), Small Parts (Canine, Feline, Equine), MSK (Canine, Feline, Equine), Vascular (Canine, Feline, Equine), Ocular (Canine, Feline), Reproduction (Canine, Feline)	NGB-053, available, multi-angle, reusable
L16-4Hs		Linear	• Bandwidth: 4.0 – 12.8 MHz • Field of View (max): 2.53 cm	Abdomen (Canine, Feline), Small Parts (Canine, Feline, Equine), MSK (Canine, Feline, Equine)	not available
L14-3Ws		Linear	• Bandwidth: 3.0 – 15.0MHz • Field of View (max): 5.08 cm	Abdomen (Canine, Feline), Small Parts (Canine, Feline),MSK (Canine, Feline), Ocular (Canine, Feline), Reproduction (Canine, Feline)	NGB-054, available , multi angle, reusable

Probe	Illustration	Type of Array	Parameters	Application	Biopsy kit
L24-6Hs		Linear	• Bandwidth: 6.0 – 24.0 MHz • Field of View (max): 2.85 cm	Abdomen (Canine, Feline), MSK (Canine, Feline), Small Parts (Canine, Feline), Ocular (Canine, Feline)	not available
P4-2s		Phased	• Bandwidth: 1.5 – 4.5 MHz • Field of View (max): 90°	Cardiology (Canine, Equine, Bovine, Ovine)	NGB-011, available, multi angle, reusable
P8-2s		Phased	• Bandwidth: 2.3 – 8.0 MHz • Field of View (max): 90°	Cardiology (Canine, Feline), Abdomen (Canine, Feline)	not available
P10-4s		Phased	• Bandwidth: 3.0 – 11.4 MHz • Field of View (max): 90°	Cardiology (Canine, Feline), Abdomen (Canine, Feline)	not available
6LE5Vs		Intrarectal	• Bandwidth: 3.0 – 11.5 MHz • Field of View (max): 4.73 cm	Reproduction (Equine, Bovine, Ovine)	not available

mindray
animal medical